LABORATORY SOLUTION

Milling & Grinding Equipments

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Main Products

FRITSCH made in Germany
The Expert of Milling to Nano-scale


www.mjr.co.kr | www.mj-research.co.kr

CONTENT

3	Sample Preparation
4	Principle of High-Energy Ball Mill
5	Line-up of High-Energy Ball Mill
6	High-Energy Ball Mill
10	Vibratory Ball Mill
12	Rotor Mill
13	Cutting Mill
16	Jaw Crusher
17	Disk Mill
18	Cup Mill
19	Hand Mortar Grinder
20	Cross Beater Mill
21	Knife Mill
22	Sieve Shaker
24	Particle Size Analyzer
25	Image Sizer
26	Sample Divider / Feeder
27	Laboratory for Sample Test

HOW to Choose the right type of comminution


Comminution principle	Type of comminution	Mill equipments
	Pressure Moving surfaces apply pressure to the sample in between.	Jaw crusher
	Impact The sample is accelerated against a surface at an extremely high speed.	Planetary ball mill, Ball mill Cup mill Rotor mill
	Friction With friction-based comminution, the sample is ground between two surfaces due to the vertical pressure of one surfaces relative to the other during simultaneous movements.	Mortar grinder, Disk mill
	Shearing Here, comminution takes place due to the shearing effect between at least on fixed surfaces and one moving surface that both move against one another.	Disk mill Rotor mill
	Cutting Rotating knives cut the sample between blades or fixed knives.	Cutting mill Knife mill

Abrasive	Planetary mill, Ball mill, P9
Alloy	Planetary mill, Ball mill, P1
Analytic	Planetary mill, Ball mill, P11, P14
Animal food	Cutting mill, Ball mill
Biology	Cutting mill, Planetary mill, Ball mill, P11
Bones	Cutting mill, P0
Building material	Planetary mill, Ball mill, P1, P13, P9
Catalytic converter	Planetary mill, Ball mill, P1
Cement clinker	Planetary mill, Ball mill, P1, P13, P9
Ceramics	Planetary mill, Ball mill, P1, P13, P9
Coal	Planetary mill, Ball mill, P1, P13
Coating powder	Planetary mill, Ball mill, P14
Compound material	Cutting mill, P14
Drug	Cutting mill, Planetary mill, Ball mill
Electronic scrap	Cutting mill, P0, P14
Environment	Cutting mill, P11, Ball mill
Feed (pellet)	Cutting mill, P11, P14
Fertilisers	Planetary mill, Ball mill, P14
Fibres	Planetary mill, Ball mill, P1, P13, P9
Foils	Cutting mill
Food	Cutting mill, Planetary mill, Ball mill, P11
Glass	Planetary mill, Ball mill, P1, P13, P9
Grains	Cutting mill, P14
Household waste	Cutting mill
Leather	Cutting mill

Lime	Planetary mill, Ball mill
Material research	Planetary mill, P7 premium, P5 premium
Mechanical alloying/activation	Planetary mill
Metallurgy	Planetary mill, P1, P9
Minerals	Planetary mill, Ball mill, P1, P13, P9
Mining	Planetary mill, Ball mill, P1, P13, P9
Ores	Planetary mill, Ball mill, P1, P13, P9
Pharmaceutical	P14
Pigments	Planetary mill, Ball mill
Plants	Cutting mill, P11, P14
Plastics	Cutting mill, P14
Refractory material	Planetary mill, Ball mill, P1, P13, P9
Rocks	Planetary mill, Ball mill, P1, P13, P9
RoHS	Cuttingm mill, P0, P9, P14
Rubber	Cutting mill, P14
Sediments	Planetary mill, Ball mill, P9
Slags	Planetary mill, P1, P9
Soil research	Planetary mill, Ball mill, P1, P13
Spectroscopy	Planetary mill, Ball mill, P14
Spices	Cutting mill, P14
Tablets	Planetary mill, Ball mill
Textiles	Cutting mill
Wood	Cutting mill, P14

* Planetary Mill: High-energy Ball Mill. Ball Mill: Vibratory Ball Mill.
** P0: Vibratory Ball Mill. P1: Jaw Crusher. P9: Cup Mill. P11: Knife Mill. P13: Disk Mill. P14: Rotor Mill.

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
Process of the Planetary Ball Mill

Line up of High-Energy Ball Mill

	Planetary Ball Mill PULVERISETTE 5 Classic line	Planetary Mono Mill PULVERISETTE 6 Classic line	Planetary Micro Mill PULVERISETTE 7 Classic line	Vario-Planetary Mill PULVERISETTE 4 Classic line	Planetary Ball Mill PULVERISETTE 5P Premium line	Planetary Micro Mill PULVERISETTE 7P Premium line
Working principle	Impact force	Impact force	Impact force	Impact force	Impact force	Impact force
Number of working stations	4	1	2	2	2	2
Grinding bowl size	80, 250, 500 ml	80, 250, 500 ml	12, 45 ml	12, 45, 80, 250, 500 ml	125, 150, 250, 420, 500 ml	20, 45, 80 ml
Grinding ball diameter	0.1 ~ 40 mm	0.1 ~ 40 mm	0.1 ~ 15 mm	0.1 ~ 40 mm	0.1 ~ 40 mm	0.1 ~ 15 mm
Max. feed size	10 mm	10 mm	5 mm	10 mm	10 mm	5 mm
Min. sample quantity	10 ml	10 ml	0.5 ml	0.5 ml	15 ml	2 ml
Max. sample quantity	900 ml	225 ml	40 ml	450 ml	450 ml	70 ml
Final fineness	< 1 um	< 1 um	< 1 um	< 1 um	< 0.1 um	< 0.1 um
Grinding process	Dry / wet	Dry / wet	Dry / wet	Dry / wet	Dry / wet	Dry / wet
Grinding in inert gas	Yes	Yes	Yes	Yes	Yes	Yes
Gas press. & temp. measurement	Yes	Yes	No	Yes	Yes	Yes
Rotation speed of main disk	50 ~ 400 rpm	100 ~ 650 rpm	100 ~ 800 rpm	0 ~ 400 rpm	100 ~ 800 rpm	100 ~ 1,100 rpm
Transmission ratio Planetary disk / grinding bowl	$i_{relative} = 1 : -2.19$	$i_{relative} = 1 : -1.82$	$i_{relative} = 1 : -2$	Variable	$i_{relative} = 1 : -2$	$i_{relative} = 1 : -2$
Effective diameter of main disk	~ 250 mm	121.6 mm	140 mm	~ 250 mm	180 mm	140 mm
Centrifugal acceleration	22 g	29 g	50 g	22 g	64 g	95 g
Interface	Yes	Yes	Yes	Yes	Yes	Yes
Electrical details	200~240V, 1730 watt	200~240V, 1000 watt	200~240V, 740 watt	200~480V, 6000 watt	200~240V, 2800 watt	200~240V, 1200 watt
Weight	120 kg	63 kg	35 kg	320 kg	110 kg	44 kg
Dimension w x d x h	580 x 670 x 570 mm	370 x 530 x 500 mm	370 x 530 x 500 mm	600 x 800 x 110 mm	820 x 520 x 480 mm	400 x 580 x 360 mm

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PULVERISETTE 7

Ideal for the smallest quantity

Working station	2 ea
Rotation speed	100 ~ 800 rpm
Max. G force	50 g
Grinding bowl	12, 45 ml
Max. sample volume	40 ml
Min. sample volume	0.5 ml
Max. feed size	5 mm
Final fineness	5 um (dry), < 1um (wet)
Inert gas	possible (optional)
Dimensions (w x d x h)	370 x 530 x 500 mm
Weight	35 kg

More safety operation

Premium line advantage

Easy insertion

ServoLOCK

Clamping release

Automatic closing

High-performance grinding down into **Nano range**

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Application

- 세라믹재료 : Porcelain, sintered ceramics, clay, etc.
- 금속재료 : Coating powder, metal, sinter, ores, etc.
- 화학 : Organic, inorganic materials, pigment, glass, etc.
- 광물 및 암석 : Rock, gravel, sand, mineral, etc.
- 제약원료 : Drug, tablet, extract, etc
- 분석시료 : XRF, Spectroscopy, XRD, etc.
- 바이오 : Plant, freeze-dried sampls, etc.
- Mechanical alloying

Ferrovanadium

Glass

Soil

Granite

Planetary Mill Premium line

PULVERISETTE 5 Premium

Ideal for large quantities

Working station	2 ea
Rotation speed	100 ~ 800 rpm
Max. G force	64 g
Grinding bowl	125, 150, 250, 420, 500 ml
Max. sample volume	450 ml
Min. sample volume	40 ml
Max. feed size	10 mm
Final fineness	5 um (dry), < 0.1 um (wet)
Inert gas	possible (optional)
Dimensions (w x d x h)	820 x 520 x 480 mm
Weight	110 kg

Planetary Micro Mill Premium line

PULVERISETTE 7 Premium

Ideal for the smallest quantity

Working station	2 ea
Rotation speed	100 ~ 1,100 rpm
Max. G force	95 g
Grinding bowl	20, 45, 80 ml
Max. sample volume	70 ml
Min. sample volume	2 ml
Max. feed size	5 mm
Final fineness	5 um (dry), < 0.1 um (wet)
Inert gas	possible (optional)
Dimensions (w x d x h)	400 x 580 x 360 mm
Weight	44 kg

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Grinding and sieving in one unit

Application

- 세라믹재료 : Porcelain, sintered ceramics, clay, etc.
- 금속재료 : Coating powder, metal, sinter, ores, etc.
- 화학 : Organic, inorganic materials, pigment, glass, etc.
- 광물 및 암석 : Soil, sand, mineral, etc.
- 제약원료 : Drug, tablet, extract, etc
- 분석시료 : XRF, Spectroscopy, XRD, etc.
- 바이오 : Hair, freeze-dried sampls, bone, Teeth, etc.
- RoHS : M.P board, camera, keypad, LCD glass panel, etc

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
Application

- 열에 민감한 샘플: PVC, PP, PE, films, plastic, etc.
- Polymer, rubber, resins, chemicals, filler, etc.
- 식품: Rice, Dried food, coffee, seed, grain, peat, etc.
- 제약: Tablet, granuals, 제약원료, etc.
- 바이오: Plants, root, leaves, etc.
- RoHS: Electronic parts, glass, etc.

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
Grinding principle	Impact force
Variable speed	6,000 ~ 20,000 rpm
Max. throughput	5 L / h
Ideal sample quantity	5 ~ 10 ml
Max. feed size	10 mm
Final fineness	40 um
Dimensions (w x d x h)	310 x 480 x 470 mm
Weight	23 kg

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Application


- Plastics & textiles: PVC, PP, PE, plastic, rubber, fabric, etc.
- Agriculture: Wood, leaves, roots, grain, corn, peat, etc.
- 환경: Paper, metal-free waste, household waste, etc.
- 제약: Tablet, granuals, 제약원료, etc.
- 바이오: Bone, plant, etc.
- 식품: Maize, malt, pasta, herbs, spices, dried meat, etc.
- RoHS: Electronic parts, glass, etc.

3-Type Rotor

Sieve cassette

Mini Cutting Mill
PULVERISETTE 29

Table-top Cutting mill


Grinding principle	Cutting
Max. feed size	up to 13 mm dia.
Sieve inserts	0.25 ~ 6 mm
Max. throughput	up to 100 ml
Feeding	Batch / continuous
Material of cutting tool	Hardened stainless steel
Rotor speed	500 ~ 6,000 rpm
Dimensions (w x d x h)	370 x 300 x 300 mm
Weight	15 kg

Cutting Mill
PULVERISETTE 15

General cutting mill

Grinding principle	Cutting
Max. feed size	40 x 70 mm
Sieve inserts	0.25 ~ 6 mm
Max. throughput	50 l / h
Feeding	Batch / continuous
Material of cutting tool	Tool steel, chromium-free
Rotor speed	2,800 , 3,400 rpm
Rotor peripheral speed	14.92 m/s, 17.9 m/s
Dimensions (w x d x h)	420 x 480 x 690 mm
Weight	42 kg

Universal Cutting Mill
PULVERISETTE 19

Fast effective comminution

Grinding principle	Cutting
Max. feed size	40 x 80 mm
Sieve inserts	0.2 ~ 6 mm
Max. throughput	60 l / h
Feeding	Batch / continuous
Material of cutting tool	Steel, T.C., chromium-free
Rotor speed	300 ~ 3,000 rpm
Rotor peripheral speed	1.57 ~ 15.7 m/s
Dimensions (w x d x h)	400 x 790 x 560 mm
Weight	60 kg

Power Cutting Mill
PULVERISETTE 19 L

Powerful pre-crushing for large sample

Grinding principle	Cutting
Max. feed size	120 x 85 mm
Sieve inserts	1 ~ 10 mm
Max. throughput	85 l / h
Feeding	Batch / continuous
Material of cutting tool	Steel, T.C., chromium-free
Rotor speed	300 ~ 3,000 rpm
Rotor peripheral speed	1.57 ~ 15.7 m/s
Dimensions (w x d x h)	450 x 650 x 630 mm
Weight	88 kg

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Grinding principle	Pressure
Max. feed size	95 mm
Final fineness	1 ~ 15 mm
Jaw size (w x d)	100 x 260 mm
Min. sample quantity	20 ml
Max. throughput	200 kg/h
Feeding	Batch / continuous
Grinding parts	Fixed & movable jaws
Eccentric oscillations	308 movements / min
Dimensions (w x d x h)	400 x 800 x 800 mm
Weight	205 kg

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Grinding principle	Shearing
Max. feed size	20 mm
Final fineness	0.1 ~ 12 mm
Disk diameter	200 mm
Min. sample quantity	20 ~ 30 ml
Max. throughput	150 kg/h
Feeding	Batch / continuous
Grinding parts	Fixed & movable disks
Rotating speed of disk	440 rpm
Dimensions (w x d x h)	440 x 870 x 400 mm
Weight	140 kg

Application

- 암석 및 토양: Clinker, quartz, concrete, etc.
- 광업: Slag, coke, coal, vanadium, ores, etc.
- Geology: Granite, silicate, other rocks, etc.
- Ceramic: Sintered ceramic, Dental ceramic, glass, etc.

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
Glass bottle 분쇄 후

Cup Mill [Disk mill, Shutter, Ring mill]

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Grinding principle	Impact
Max. feed size	12 mm
Final fineness	10 ~ 20 um
Min. sample quantity	10 ~ 20 ml
Max. sample quantity	250 ml
Feeding	Batch
Grinding parts	Grinding puck with ring
Motor speed	1,100 rpm, 1,300 rpm
Useful capacity	50, 100, 250 ml
Dimensions (w x d x h)	720 x 660 x 1,200 mm
Weight	250 kg

Application

- Mining: Coal, ores, minerals, etc.
- Metallurgy: Slag, cast iron sample, etc.
- Agriculture & environment: Soil, sludge, etc.
- Analysis: XRF, etc.

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
Application

- Agriculture and forestry: Grain, peat, seeds, dried plants
- Geology: Salts, gypsum, potash, stons, soil, limestone, etc
- Mining: Coal, coke, ores, minerals, slate, slag, etc.
- Plastics and textiles: Resins, cellulose, etc.
- Ceramics: Oxide ceramics

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
Grinding principle	Impact
Max. feed size	20 mm
Min. sample quantity	30 ~ 40 ml
Max. throughput	80 L / h
Sieve insert	0.12 ~ 10 mm
Feeding	Batchwise / continuous
Grinding part	Cross beater
Materials of grinding parts	Cast iron, Stainless steel
Motor speed	2,850 rpm
Dimensions (w x d x h)	420 x 450 x 560 mm
Weight	38 kg

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Application

- 식품: Fish, meat, ham, vegetable, potato, cheese, grain seed, nut, chocolate, bread, etc.
- Animal food: Feed pellet, dog food, etc.
- 제약: Tablet, granuals, 제약원료, etc.
- 바이오: Plant, leaves, freeze-dried samples, etc

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
Vibratory Sieve Shaker

ANALYSETTE 3 PRO

Precise sieving with automatic amplitude control

Sieving principle	Vibrating
Dry sieving	20 um ~ 63 mm
Wet sieving	20 um ~ 10 mm
Micro-sieving	5 um ~ 100 um
Max. weight of sieve stack	3 kg
Amplitude	0.1 ~ 3 mm
Amplitude control	automatic
Sieve diameter	75 mm, 100 mm, 203 mm
Max. number of sieve	10 (50 mm ht.), 16 (25 mm ht.)
Dimensions (w x d x h)	370 x 400 x 200 mm
Weight	21 kg

Heavy Duty Analytical Sieve Shaker

ANALYSETTE 18

Effective sieving of large quantity

Sieving principle	Vibrating
Dry sieving	20 um ~ 125 mm
Max.weight of sieve stack	42 kg
Amplitude	0.1 ~ 2 mm
Amplitude control	Automatic
Sieve diameter	200, 250, 300, 315, 350, 400, 450 mm
Max. number of sieve	12 (65 mm height)
Max. height of sieve stack	845 mm
Dimensions (w x d x h)	580 x 590 x 1,300 mm
Weight	135 kg

Particle Size Analyzer

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
Feature


- 레이저회절법을 이용한 입도측정
- 입도측정범위: Micro 0.5 ~ 1,500 um
Nano 0.01 ~ 3,800 um
- Fraunhofer diffraction or Mie scattering
- Green laser를 이용하여 입도분석
- Wider scattering angle range
- Wet: Sample vol. 150~500 ml
Ultrasonic suspension 50 wat
- Automatic dispersion, measurement
Analysis, rinsing, reporting
- 구성: Measuring unit, Wet dispersion unit,
Ultrasonic box
- All method SOP control

Principle



Laser diffraction method (레이저회절법)

Measurement design
ANALYSETTE 22 NeXT Micro

Measurement design
ANALYSETTE 22 NeXT Nano

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Image Analysis of Particle shape and size

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Measuring range	Dry : 20 um ~ 20 mm, Wet : 5 um ~ 3 mm
Method of analysis	Dynamic image analysis
Type of analysis	Dry : powder, solid Wet : suspension, emulsion
Measurement values	Particle shape and size
Standard	ISO 13322-2
Camera	5 megapixel CMOS camera 2,448 x 2,050 pixel resolution
Sample quantity	Dry: 10 ~ 100 g, Wet: 0.1 ~ 1 g
Measuring speed	Max. 75 images / s
Dimensions (w x d x h)	900 x 300 x 550 mm
Weight	Dry: 36.8 kg, Wet: 58.8 kg
System	Window10, 16GB RAM

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About FRITSCH

1920

Company founded as a technical precious stones dealer

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MJ Research Co., Ltd.